

For more information



Visit www.advansta.com/products/HRP_conjugates or go directly to the web page by scanning the QR-code with your favorite bar-code scanner app on your smart phone.

HRP-conjugated Secondary Antibodies

HRP-conjugated secondary antibodies for chemiluminescent immunodetection

Description

Advansta's high-quality HRP-conjugated secondary antibodies offer unparalleled sensitivity and performance for immunoblotting and ELISA applications. These secondary HRP conjugated antibodies are optimized with Advansta's protein detection systems including WesternBright® ECL, WesternBright® Quantum®, WesternBright® Sirius® and ELISABright™.

Advansta Corporation

2140 Bering Drive | San Jose, CA 95131

Tel: 650.325.1980 | Fax: 650.325.1904

Email: support@advansta.com

www.advansta.com



April 16, 2020 | D-26033-403

HRP-conjugated Secondary Antibodies

Advantages

- **Performance** – excellent signal in immunoblotting and ELISA applications
- **Convenience** – each antibody is demonstrated to provide optimal results with Advansta's protein detection systems
- **Purity** – each antibody is affinity purified to maximize specificity

Tips

- Western Blotting and ELISA require optimization of both primary and secondary antibody concentrations. These must be determined for every target-antibody pair.
- A good starting dilution of the HRP-conjugated secondary antibodies is 1:20,000 when using Advansta's protein detection reagents.

Ordering Information

Catalog Number	Product	Size
R-05071-500	Goat anti-mouse IgG (H+L), HRP conjugate	500 µl
R-05072-500	Goat anti-rabbit IgG (H+L), HRP conjugate	500 µl
R-05073-500	Goat anti-human IgG (H+L), HRP conjugate	500 µl
R-05074-500	Goat anti-chicken IgY(IgG) (H+L), HRP conjugate	500 µl
R-05075-500	Goat anti-rat IgG (H+L), HRP conjugate	500 µl
R-05076-500	Goat anti-guinea pig IgG (H+L), HRP conjugate	500 µl
R-05077-500	Donkey anti-goat IgG (H+L), HRP conjugate	500 µl

